



## HisExpress™ Purification Protocol for Native His-tagged Proteins from *E. coli* Cells (For research use only)

### HisExpress Kit components



### Simplified Protocol



- ① Using the syringe and tubing tip, pre-equilibrate the OmniLyse® cartridge and HisExpress™ column with 1 mL of Binding Buffer each. Take a cell pellet (from up to 20 mL of *E. coli* culture) and resuspend in 1 mL of Binding Buffer. Attach tubing tip and syringe to OmniLyse® cartridge. Plug the OmniLyse® cartridge blue plug into socket on battery pack. Draw the sample into cartridge displacing air. Turn on the battery pack and draw the remaining sample through the cartridge. Reverse direction, dispensing the sample back through the cartridge into the original sample container. Continue drawing and dispensing the sample for 3 minutes, then, turn off the battery pack and withdraw the entire sample into syringe.
- ② Attach the syringe to the HisExpress™ column and dispense the lysate through the column collecting the lysate flow-through in a 1.5 mL microfuge tube. Save for analysis or discard.
- ③ Attach a syringe containing 3 mL of Wash Buffer and dispense it through the column. Save for analysis or discard.
- ④ Attach syringe containing 1 mL of Elution Buffer and dispense through the column collecting desired sized fractions.

Product use restrictions:

OmniLyse® cartridges and HisExpress™ columns were developed and are sold for research purposes only and are suitable for *in vitro* applications only.

Limited Use Label License: Claremont BioSolutions LLC

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control.

Claremont BioSolutions LLC will not assert a claim against the buyer of infringement of intellectual property owned or controlled by Claremont BioSolutions LLC which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Claremont BioSolutions LLC is willing to accept return of the product with a full refund. For information about purchasing a license to use this product or the technology embedded in it for any use other than for research use please contact Claremont BioSolutions LLC, 1182 Monte Vista #11, Upland, California 91786, 909-920-3887.